

WidowX AI

Leader

Parts List


Parts:

- (1) *WidowX AI Leader Robot Arm*
- (1) *WidowX AI Controller Box*
- (1) *Aviator to XT30(2+2) Cable*
- (1) *24V 25A Power Supply*
- (1) *M12 X-Code 8-Pin Male to RJ45 Male Cable 1m*
- (1) *Base Cable Lock*
- (1) *Left Leader Handle*
- (1) *Right Leader Handle*
- (1) *Leader Right Index Finger Paddle*
- (1) *Leader Right Thumb Finger Paddle*
- (1) *Leader Left Index Finger Paddle*
- (1) *Leader Left Thumb Finger Paddle*
- (1) *M2 Hex Key*
- (1) *M2.5 Hex Key*
- (1) *M3 Hex Key*
- (1) *2mL Loctite Tube*
- (1) *Base Clamp Plate*
- (12) *M3x6mm Socket Head Bolt*
- (2) *M3x6mm Button Head Bolt*
- (4) *M3x8mm Button Head Bolt*
- (4) *M5x8mm Button Head Bolt*

1) Attach Arm to Base Clamp Plate

Parts:

- (1) WidowX AI Follower Robot Arm
- (1) Base Clamp Plate
- (4) M5x8mm Button Head Bolts
- (1) M3 Hex Key
- (1) 2mL Loctite Tube




2R) Attach Finger Paddles to Arm (Right Handed)

For a right-handed arm, Attach the following parts.

Parts:

- (1) WidowX AI Follower Robot Arm
- (1) Leader Right Index Finger Paddle
- (1) Leader Right Thumb Finger Paddle
- (12) M3x6mm Socket Head Bolt
- (1) M2.5 Hex Key
- (1) 2mL Loctite Tube



*Mind Orientation of
Finger Paddles-
Letters should not
be visible when mounted*



3R) Attach Handle to Arm (Right Handed)

Parts:

- (1) WidowX AI Follower Robot Arm
- (1) Right Leader Handle
- (4) M3x8mm Button Head Bolt
- (1) M2 Hex Key
- (1) 2mL Loctite Tube - Mind Orientation of Handle



2L) Attach Finger Paddles to Arm (Left Handed)

For a left-handed arm, Attach the following parts.

Parts:

- (1) WidowX AI Follower Robot Arm
- (1) Leader Left Index Finger Paddle
- (1) Leader Left Thumb Finger Paddle
- (12) M3x6mm Socket Head Bolt
- (1) M2.5 Hex Key
- (1) 2mL Loctite Tube


*Mind Orientation of
Finger Paddles-
Letters should not
be visible when mounted*



3L) Attach Handle to Arm (Left Handed)

Parts:

- (1) WidowX AI Follower Robot Arm
- (1) Left Leader Handle
- (4) M3x8mm Button Head Bolt
- (1) M2 Hex Key
- (1) 2mL Loctite Tube - Mind Orientation of Handle



4) Connect Cables to Controller Box

Parts:


- (1) WidowX AI Controller Box
- (1) Aviator to XT30(2+2) Cable
- (1) 24V 25A Power Supply
- (1) M12 X-Code 8-Pin Male to RJ45 Male Cable 1m



5) Connect Cable to Arm

Parts:

- (1) WidowX AI Follower Robot Arm
- (1) Aviator to XT30(2+2) Cable
- (1) Base Cable Lock
- (2) 3x6mm Button Head Bolt
- (1) M2 Hex Key
- (1) 2mL Loctite Tube



6) Power on Controller Box

Parts:

- *WidowX AI Follower Robot Arm*
- *WidowX AI Controller Box*

With Cables attached, turn on controller and wait for the light to turn green. You are now ready to use your WidowX AI Robot Arm.

